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
参考文献：1. World Health Organization (WHO). Naming the coronavirus disease (COVID-19) and the virus that causes it. データアクセス日：2020年3月30日。Accessible here: [https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-\(covid-2019\)-and-the-virus-that-causes-it](https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance/naming-the-coronavirus-disease-(covid-2019)-and-the-virus-that-causes-it) 2. World Health Organization (WHO). Clinical management of severe acute respiratory infection (SARI) when COVID-19 disease is suspected: Interim guidance V 1.2. データアクセス日：2020年3月30日。Accessible here: <https://apps.who.int/iris/bitstream/handle/10665/331446/WHO-2019-nCoV-clinical-2020.4-eng.pdf?sequence=1&isAllowed=y> 3. Peng QY, Chinese Critical Care Ultrasound Study Group (CCUSG), et al. Findings of lung ultrasonography of novel coronavirus pneumonia during the 2019–2020 epidemic. Intensive Care Med. Published online 12 March 2020. DOI : <https://doi.org/10.1007/s00134-020-05996-6> 4. Buonsenso D et al. Lancet Respir Med. Published online 20 March 2020. DOI: [https://doi.org/10.1016/S2213-2600\(20\)30120-X](https://doi.org/10.1016/S2213-2600(20)30120-X) 5. Spaulding EH. Chemical disinfection of medical and surgical materials. In: Lawrence C, Block SS, editor. Disinfection, sterilization, and preservation. Philadelphia (PA): Lea & Febiger; 1968. p. 517-31. 6. Centre for Disease Control and Prevention (CDC) 2008. Guideline for Disinfection and Sterilization in Healthcare Facilities. 7. Costello C et al. Prevention of pathogen transmission during ultrasound use in the Intensive Care Unit: Recommendations from the College of Intensive Care Medicine Ultrasound Special Interest Group (USIG). Published online 26 March 2020. DOI: <https://doi.org/10.1002/ajum.12205> 8. The International Society of Ultrasound in Obstetrics and Gynecology (ISUOG) Safety Committee Position Statement: safe performance of obstetric and gynecological scans and equipment cleaning in the context of COVID-19. データアクセス日：2020年3月30日。Accessible here: <https://www.isuog.org/resource/isuog-safety-committee-position-statement-safe-performance-of-obstetric-and-gynecological-scans-and-equipment-cleaning-in-the-context-of-covid-19.html> 9. Wessner CE et al. Journal of Diagnostic Medical Sonography. September 2020. doi:10.1177/8756479320959035 10. Sheth S et al. Guidelines for Ultrasound in the Radiology Department During the COVID-19 Pandemic. Ultrasound Quarterly. September 2020;36(3):200-205.

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